

MALANG

MalangMassageChair.com

Customer support



Scan for video guide
and latest manual

© 2026 MALANG. All Rights Reserved.

MALANG

Massage Chair

MODEL ML-FV8-25V1

User Manual



————— Designed for Daily Recovery —————

12. Troubleshooting Guide

This section provides guidance for identifying and resolving common operating conditions associated with your massage chair.

Operating Sounds (Normal)

Possible Causes	Solution
Sound from the massage mechanism during operation	These sounds are a normal part of operation and do not indicate a malfunction. No action is required.
Airbags inflating and deflating	
Solenoid valve activation	
Air pump operation	

The Chair Does Not Operate

Possible Causes	Solution
The power plug is not fully inserted	Ensure the power plug is securely connected to both the wall outlet and the chair
The main power switch is turned off	Confirm that the main power switch is turned on
The power button has not been activated	Press the power button on the remote control or shortcut panel

Massage Stops During Use

Possible Causes

- The 20-minute session limit has been reached and the chair has entered standby mode
- The motor has overheated and activated protective shutdown

Solution

- Press the power button to restart the massage session
- Turn off the chair and allow the motor to cool before resuming use

Important Notes

- The manufacturer reserves the right to improve the product and this manual without prior notice.
- Product appearance may vary slightly from the images shown in this manual.
- While this manual has been carefully reviewed, minor printing or interpretation errors may occur.
- Please dispose of this product in accordance with local environmental regulations.

Welcome

Welcome to refined comfort.

This massage chair is precisely engineered to support the natural movement and balance of the body. Each contour, motion, and function is designed to deliver a seamless full-body experience—quiet, fluid, and deeply restorative.

Gentle warmth releases tension. Air compression moves in slow, rhythmic waves. Subtle vibration aligns naturally with your breathing, guiding the body toward a state of calm and renewal.

Every session is a moment of intentional rest—a private ritual created to restore comfort, clarity, and balance.

Please read all safety instructions carefully before use to ensure optimal performance and long-term reliability.

Keep this manual for future reference.

Contents

01. Product Features
02. Component Names
03. How to Use
04. Tablet Interface
05. Technical Specifications
06. Maintenance
07. Safety Precautions
08. Safety Protections
09. Important Tips
10. Product Warranty
11. After-Sales Service & Warranty Card
12. Troubleshooting Guide

Massage Duration

The recommended maximum session time is 20 minutes. For safety, the chair will automatically stop after extended continuous use.

Heat Function

Use the heat function with caution, especially if you are sensitive to warmth.

Sensitive Areas

Do not apply massage to sensitive areas such as elbows, knees, chest, or abdomen. If discomfort occurs, stop use immediately and consult a medical professional.

Power Interruptions

In the event of a power outage, exit the chair carefully. Manual repositioning may be required.

10. Product Warranty

To protect consumer interests, this product is covered by a two-years limited warranty, in accordance with applicable national regulations.

10.1 Complimentary Warranty Service

If a product malfunction occurs due to a manufacturing or material defect within the two-year warranty period, necessary replacement parts will be provided at no charge.

10.2 Warranty Exclusions (Chargeable Service)

The warranty does not cover damage or malfunction resulting from the following conditions. In such cases, service and replacement parts may be provided at the customer's expense:

- Improper use, misuse, or unauthorized self-repair
- Use of a power supply that does not meet specified voltage requirements
- Damage caused by unauthorized relocation, transportation, modification, or installation
- Damage resulting from normal wear, contamination, or lack of routine care
- Damage caused by natural disasters, environmental events, or pest infestation

10.3 Out-of-Warranty Service

Repair services will be provided at full cost under the following circumstances:

- The warranty period has expired
- A valid warranty card cannot be presented or registered


SL-Type Extended Guide Rail

The advanced SL-track allows the dual massage mechanisms to follow the body's natural curvature, delivering synchronized, full-length coverage from the neck through the thighs for a continuous and balanced massage experience.

Immersive Hi-Fi Bluetooth® Audio

Integrated high-fidelity speakers deliver clear, enveloping sound, enhancing relaxation and creating a refined, immersive sensory environment.

2. Component Names



3. How to Use

3.1 Installing the Power Cord

Step 1

Locate the power inlet on the rear of the massage chair (Figure 1).

Step 2

Firmly insert the power cord plug into the power inlet.

Step 3

Connect the opposite end of the power cord to a standard electrical outlet.

3.2 Moving the Massage Chair

Step 1

Ensure the chair is unoccupied and powered off.

Step 2


Stand behind the chair and press down firmly on the backrest.

Step 3

Tilt the chair forward until it rests securely on the moving wheels (approximately 40°).

Step 4

Carefully push the chair forward to its new position.



Important

- Do not apply excessive force when tilting the chair, as this may cause instability.
- Never move the chair while it is occupied.
- To protect flooring, place a mat or protective padding beneath the chair, particularly during relocation.

Intended Use

This massage chair is designed for personal, in-home use only. It is not intended for commercial, rental, or medical applications.

Temperature

If the chair has been stored in a low-temperature environment, allow it to reach room temperature for at least one hour before use to ensure proper operation.

9. Important Tips

Jewelry and Clothing

Remove jewelry such as rings, bracelets, necklaces, and pins before use. Avoid loose or restrictive clothing that may become caught during operation.

Initial Use

During first use, only the power function is available. Additional functions will activate once the chair has fully powered on.

Children and Supervision

Children should use the massage chair only under adult supervision. Do not allow children to climb, play, or stand on the chair.

Massage Intensity

Begin at a low intensity and increase gradually according to comfort. Massage pressure should feel relaxing and supportive, never painful.

Health Considerations

Consult a qualified medical professional before using this product if you have a medical condition. Do not use the chair if you have open wounds, skin conditions, heart disease, high blood pressure, or other serious health concerns without medical guidance.

Electrical Safety

Do not use the chair in combination with other therapeutic electrical devices, such as electric blankets. Avoid use within one hour after consuming a full meal.

Pinch Points

Keep hands and body parts away from moving components and areas where compression or pinching may occur.

Flammable Materials

Keep all flammable or combustible materials away from the chair at all times.

Supervision

When used by children, elderly individuals, or persons with physical or cognitive limitations, close adult supervision is required.

Authorized Parts

Use only replacement parts and accessories approved by the manufacturer.

Damage or Malfunction

Do not operate the chair if the power cord or plug is damaged, if the chair has been dropped, exposed to water, or is not functioning properly. Contact an authorized service center for inspection and repair.

Moving the Chair

Do not drag the chair by the power cord or use the cord as a handle.

Heat Exposure

Keep the power cord away from heating surfaces.

Air Ventilation

Ensure that air outlets remain unobstructed and free from dust, hair, or debris.

Foreign Objects

Do not insert any objects into the openings of the chair.

Location

Use the chair indoors only. Do not operate in outdoor environments or in areas where aerosol products or oxygen are in use.

Power-Off Procedure

After turning off all operating controls, unplug the power cord from the wall outlet.

Coverings

Do not cover the massage chair during use or operation.


8. Safety Protection

Environment

Operate the massage chair in a clean, dry, and well-ventilated space. Avoid excessive humidity, dust, direct sunlight, or proximity to heat sources, as these conditions may affect performance or electrical safety.

4. Tablet Interface

The tablet interface provides intuitive control over all massage functions. Its streamlined layout allows effortless customization while maintaining a calm and immersive user experience.



4.1 Body Scanning

The chair automatically scans the body to detect key pressure points and alignment, allowing the massage system to adapt precisely to individual body contours for a personalized session.

Heating Therapy

Activate or deactivate heat therapy for the waist and leg areas to enhance muscle relaxation and circulation.

Foot Roller Massage

Adjust the intensity of the foot rollers to suit personal comfort and preference.

Leg Kneading

Select and customize leg kneading modes to deliver targeted, rhythmic lower-body relief.

Additional Controls

Preferences

Access the settings menu to adjust the following options:

- Massage timer
- Screen brightness
- Button and media volume
- Bluetooth® on/off
- Ambient lighting control
- Language selection
- Voice control activation

Pause / Resume


Press to pause an active massage session. Press again to resume.

Gentle Awakening (Stretching)

Initiates a full-body stretching program designed to gently awaken and extend the body.

More Modes

Access additional massage programs to further personalize your experience.



6. Maintenance

Cleaning

For daily care, wipe the chair with a soft, dry cloth. If necessary, use a mild, neutral detergent. Avoid corrosive, abrasive, or solvent-based cleaning agents, as they may damage the surface.

Power Cord

If the power cord becomes damaged, discontinue use immediately and contact an authorized service provider for replacement.

Electrical Safety

To reduce the risk of accidents, unplug the chair when not in use and before cleaning or servicing.

Protection

Do not place sharp or heavy objects on the chair, as this may cause permanent damage.

Storage

Store the chair in a dry, well-ventilated environment. Excessive humidity may damage internal components.

Handling

Avoid applying excessive force or pressure to any part of the chair during use or relocation.

7. Safety Precautions

Please read and follow the instructions below to ensure safe operation and long-term performance.

Power Supply

Use only an AC 110V power source.

Unplugging

Always unplug the chair from the power outlet when not in use, before cleaning, or before installing or removing components.

Overheating Prevention

Do not place cushions, blankets, or coverings on the chair during operation, as this may restrict ventilation and cause overheating.

4.8 Shortcut Instructions

The shortcut control panel, located on the side of the chair, provides quick and convenient access to frequently used functions.

Shortcut Functions

Power

Press to turn the massage chair on or off. Rotate the control to adjust the massage speed.

Auto

Press to cycle through the available automatic massage programs.

Heating

Activate or deactivate heat therapy for the waist and leg areas.

Airbag

Control the air massage system. Press to activate sequential inflation and deflation for the shoulders, arms, legs, and feet.

Roller

Turn the foot roller massage on or off.

Backrest Up

Press and hold to raise the backrest to an upright position.

Backrest Down

Press and hold to recline the backrest to a comfortable position.

5. Technical Specifications

Specification	Details
Product Name	MALANG Massage Chair ML-FV8-25V1
Product Type	Dual Roller SL-Track Massage Chair
Rated Voltage	AC 110V
Frequency	60 Hz
Rated Power	150W

4.2 Automatic Massage Programs

Indulge in a curated collection of expertly designed automatic massage programs. Each mode is thoughtfully crafted to support balance, recovery, and deep relaxation—allowing body and mind to return to harmony with a single touch.

(Program names displayed on the screen may vary by model and software version.)

Immerse yourself in pure tranquility with Serenity Mode or drift effortlessly into deep rest through Dream Flow. Experience precise relief with Upper Relief or the rejuvenating alignment of Chiropractic Care.

Awaken your senses with Gentle Awakening (Stretching) or revitalize fatigued muscles in Active Recovery. Reconnect your body's natural rhythm through Balance Yoga, embrace traditional healing with Tuina Therapy, or feel the soft embrace of air pressure in Air Embrace.

For everyday comfort, unwind with Focus Relax or enjoy soothing ease in Easy Watching, a perfect companion for your favorite moments of rest.

Each program is thoughtfully designed to deliver a sensory journey of balance, renewal, and deep relaxation—transforming every session into a refined wellness ritual.

4.3 Manual Mode

Manual Mode provides precise, hands-on control, allowing you to tailor each massage to your individual preferences.

Massage Techniques

Select from six professional massage techniques:

- Kneading
- Rhythmic Tapping
- Kneading & Tapping
- Percussion
- Zen Shiatsu
- Bi-Directional Kneading

Massage Position

Adjust the massage focus to suit your needs:

- **Focus Massage:**
Concentrates massage action on a single targeted area.
- **Partial Relax:**
Targets a specific region between the upper and lower back.

• Upper & Lower Body

Control the upper and lower massage mechanisms independently. Selecting the upper mechanism focuses massage on the upper body, while selecting the lower

Advanced Adjustments

For further personalization, the following settings may be adjusted:

• Massage Width

Select from Wide, Middle, or Narrow massage paths. Available for flapping and tapping techniques.


• Mechanism Speed

Adjust the operating speed of both upper and lower massage mechanisms.

• Air Intensity

Control the intensity of air compression during the massage session.

With these options, each session can be precisely customized—creating a massage experience aligned with your body's unique needs and preferences.



4.4 Health Scan

The Health Scan feature provides general wellness indicators by measuring heart rate, blood oxygen level, and microcirculation.

4.7 Settings Menu

The Settings menu allows you to personalize the massage chair experience with precision and ease.

Display & Sound

Screen Brightness

Adjust the brightness of the tablet display.

Button Sound

Control the volume of interface button sounds.

Media Volume

Adjust the audio level of music played through the Bluetooth® speakers.

Connectivity & Atmosphere

Bluetooth®

Enable or disable Bluetooth® audio connection.

Ambient Lighting

Turn the chair's integrated ambient lighting on or off to enhance the surrounding atmosphere.

Voice Control

Activate or deactivate voice command functionality.

Session & Language

Timer

Set the massage duration in five-minute increments, from 5 to 30 minutes.

Language

Select the user interface language. Available options: English, Korean, Spanish.

Voice Control Notice

If no voice command is detected within 10 seconds after activation, the voice assistant will automatically enter standby mode. To resume voice control, repeat one of the activation phrases.

Heat on	OK, heating on
Heat off	OK, heating off
Zero Gravity on	OK, Zero Gravity on
Zero Gravity off	OK, Zero Gravity off
Increase strength	OK, stronger massage
Decrease strength	OK, softer massage
Faster speed	OK, rollers faster
Slower speed	OK, rollers slower
Add time	OK, extended 10 minutes
Volume up	OK, volume up
Volume down	OK, volume down
Help	Here are some commands you can use... (change the tablet screen to voice control page)
Check blood pressure	Measuring blood pressure
Lights on	Lights on
Lights off	Lights off
Good Morning	Good morning! Let's make today amazing and full of positive energy.
Good Night	Good night! Wishing you sweet dreams, Rest easy and wake up refreshed!



How to Perform a Health Scan:

- Access the Health Scan interface from the tablet control.
- Place your fingertip gently on the designated glass sensor.
- Ensure full contact and remain still to avoid interference during measurement.
- Select Start Collection and wait briefly while the scan completes.
- Results for heart rate (BPM), blood oxygen level, and microcirculation will be displayed on the screen.

Note

The information provided by this feature is intended for reference only and is not designed for medical diagnosis or treatment.

4.5 Adjusting Sitting Position

The massage chair allows precise adjustment of seating position to enhance comfort and support during use.

Backrest Adjustment

Press and hold the control to raise or recline the backrest to your preferred position.

Leg Rest Extension


Press and hold to extend or retract the leg rest, allowing adjustment to suit your height.

Leg Rest Lift

Press and hold to raise or lower the entire leg rest for improved lower-body alignment.

Zero-Gravity Positions

Select from three zero-gravity reclining levels (G1, G2, G3) to reduce spinal pressure and promote deeper relaxation.



4.6 Smart Voice Control

Your massage chair is equipped with an integrated voice assistant for hands-free operation.

Activating Voice Control

Activate the voice assistant by saying one of the following phrases:

- “Hey Chair” “Hello Chair” “Hello Massage Chair”

Once activated, voice commands may be used to control massage functions, modes, and settings.

Voice Command Timeout

If no command is detected within 10 seconds after activation, the voice assistant will enter standby mode. To reactivate, repeat one of the activation phrases.

Supported Voice Commands

Voice control supports a range of functions, including:

- Power on/off
- Switching massage modes and techniques
- Adjusting massage intensity and speed
- Activating specific programs such as Serenity Mode, Deep Calm, Dream Flow, Balance Yoga, Gentle Awakening, Power Deep, and Energy Reboot
- Controlling airbag and foot roller functions
- Adjusting volume, lighting, and session duration

Turn on	OK, power on
Turn off	OK, shutting down
Next mode	OK, switching mode
Next technique	OK, switching technique
Serenity Mode	Starting Serenity Mode
Deep Calm	Starting Deep Calm
Dream Flow	Starting Dream Flow
Balance Yoga	Starting Balance Yoga
Gentle Awakening (Stretching)	Starting Gentle Awakening(Stretching)
Power Deep	Starting Power Deep
Energy Reboot	Starting Energy Reboot
Active Recovery	Starting Active Recovery
Dynamic Pulse	Starting Dynamic Pulse
Core Balance	Starting Core Balance
Chiropractic Care	Starting Chiropractic Care
Zen Touch	Starting Zen Touch
Airbags on	OK, airbags on
Airbags off	OK, airbags off
Foot roller on	OK, foot roller on
Foot roller off	OK, foot roller off
Back up	OK, raising backrest
Back down	OK, lowering backrest
Legs up	OK, raising calf
Legs down	OK, lowering calf